

### **13. Templates for interpretation of disc diffusion tests**

Zone diameters in routine tests are most easily interpreted by use of templates of zone diameter breakpoints. A program to prepare customised templates has been written and supported by Trevor Winstanley on behalf of the BSAC Working Party on Antimicrobial Susceptibility Testing.

The BSAC disc diffusion template programme is available on the BSAC web site at [www.bsac.org.uk/susceptibility-testing/bsac-disc-diffusion-template-program.cfm](http://www.bsac.org.uk/susceptibility-testing/bsac-disc-diffusion-template-program.cfm)